Heart Rate Plus by PVDApps

- a) App measures Heart rate and keeps the past history.
- b) App by PVDapps startup.
- c) Instantaneous and continuous measurement.

d)

- High accuracy
- Low error

e)

- High latency
- Fails sometimes.

Health Watcher

a) App measures Heart rate, Blood pressure, Respiration Rate, Oxygen Saturation Level.

b)

- A. . K. Kanva, C. J. Sharma and S. Deb, "Determination of SpO2 and Heart-rate using Smartphones, International Conference on Control, Instrumentation, Energy & Communication, New Delhi, India, 2014.
- extracting heart rate and respiration rate using a cell phone camera
- measuring vital signs using smart phones by b. vikram chandrasekaran.
- c) Requires less computation power and by using image processing algorithms by intel.

d)

- Uses the same data for measuring HR, BP, SpO2.
- Can be run on low end device

e)

- Fails to multiple times before giving results
- Has a high latency
- Error Prone to ambient light
- High error variance.

Instant Heart Rate

- a) App measures Heart rate very precisely and further uses data to give more insights about stress levels.
- b) **Azumio** is a mobile health company that specializes in biometric mobile technology. Azumio released Instant Heart Rate in October 2010
- c) This app has been used for research at Heart Research at UCSF(University of California San Francisco) for cardiology research training.

d)

- High accuracy
- Can be run on low end device
- Low latency
- Low run time
- Low error


e)

• Error is mostly because of instability of hand.

Plots:


Time series Data Approach with 3 seconds

- a) As we increase the sliding window the accuracy of the estimator increases.
- b) For small sliding windows error increases.



FFT

- a) Also depends upon sliding window
- b) As we increase the sliding window accuracy of the estimator increases.



TASK 4:

- 1. FPS can not less than fs = 2*f_high of the signal. f_high = 1.5. In our case i.e fs = 2*1.5 = 3 fps. As per Nyquist sampling theorem i.e. FPS > 3
- 2. The parts directly under the flash light don't change much also at the edges there are very small changes. We can choose a window at the center.
- 3. If we lower resolution we have to sample at high fps. And we sample at low fps we have to have high resolution to have meaningful results.